
ENGLISH	1 - 8
DEUTSCH	9 - 18
FRANÇAIS	19 - 26
ESPAÑOL	27 - 36
ITALIANO	37 - 44
PORTUGUÊS	45 - 52
РУССКИЙ	53 - 62


CE

1

2

**ENGLISH**

- A - Aspiration button – PP
- B - Dispense button – PP
- C - Suction speed switch – PP
- D - Dispense mode switch – PP
- E - Display
- F - Nosepiece – PP
- G - Pipette holder – SI
- H - Membranefilter – PP / PTFE
- K - Connector gasket
- M - Charging stand

NiMH battery  
Casing – PP

**FRANÇAIS**

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



6

**TABLE OF CONTENTS**

1 - INTRODUCTION .....	2
2 - WORK SAFETY INSTRUCTIONS .....	2
3 - LIMITATIONS OF USE .....	3
4 - TO SWITCH ON THE ASSISTBOY .....	3
5 - ASPIRATING AND DISPENSING LIQUIDS .....	3
6 - TROUBLESHOOTING .....	4
7 - REPLACING THE FILTER .....	6
8 - CHARGING THE BATTERY .....	6
9 - MAINTENANCE .....	7
10 - COMPONENTS .....	7
11 - ORDERING INFORMATION .....	8
12 - SPARE PARTS .....	8
13 - PRODUCT DISPOSAL .....	8

**1 - INTRODUCTION**

The **Assistboy** pipette controller is a device intended for pipetting liquids with the use of measuring pipettes. It can work with all types of glass or plastic pipettes in the volume range from 0.5 ml to 100 ml.

Two dispense modes permit selection of dispensing intensity depending on the user's needs. The selected setting of the pipette controller mode is shown on the display (fig. 1E)

Fig. 1 shows the external parts of the pipette controller with a description of the materials used.

Display indicators are shown in fig. 2.

**2 - WORK SAFETY INSTRUCTIONS**

Before starting the work with the **Assistboy** every user should read these operating instructions carefully.

Using the device inconsistently with the operating instructions may result in damaging the device.

The device should be maintained only at an authorised service centre, otherwise the manufacturer will be relieved from any liability under the warranty.

- During the work with the **Assistboy** general safety regulations regarding risks related with laboratory work should be observed. Protective clothing, goggles and gloves should be worn.
- The **Assistboy** shall be used only for measuring liquids in conditions specified by the manufacturer, which are limited due to the chemical and mechanical resistance of the device, as well as the user safety.
- The pipette controller should not be used in an environment where explosion risk is present.
- The information and instructions provided by the manufacturers of the reagents must be observed.
- Only original spare parts and accessories, recommended by the manufacturer, shall be used.
- Only the original charger, supplied by the manufacturer, shall be used for charging the battery.
- In case of incorrect functioning of the pipette controller, work shall be stopped. The device shall be cleaned in accordance with the operating instructions and/or sent for repair to an authorised service centre.
- In the case of mechanical damage to the casing, the device shall be immediately sent for repair to an authorised service centre.
- The use of excessive force during work shall be avoided.

### 3 - LIMITATIONS OF USE

- The **Assistboy** may not be used for measuring substances, the vapours of which damage the following plastics: PP, SI, ABS, EPDM, POM.
- The pipette controller may not be used in an environment where explosion risk is present.
- Flammable liquids shall not be measured - in particular substances with flash-point below 0°C (ether, acetone).
- The pipette controller shall not be used for drawing acids with a concentration above 1 mol.
- The pipette controller shall not be used for drawing solutions with a temperature above 50°C.
- The pipette controller may work in temperature range from +10°C to +35°C.

The pipet controller is suitable for general laboratory use only. It must be used only by personnel who know the health risks associated with the substances that are normally used with this instrument.

### 4 - TO SWITCH ON THE ASSISTBOY

The **Assistboy** pipette controller is switched on by pressing the aspiration button (fig. 1A) or the dispense button (fig. 1B). The display will show the selected aspiration mode, dispense mode and the battery level indicator. Examples of the display indications are shown in (fig. 5). The battery is discharged and requires charging if the indicator displays one "bar" (fig. 5.3). When the battery is fully charged, the indicator displays three „bars" (fig. 5.1). The **Assistboy** switches off automatically if not used for 5 minutes.

The **Assistboy** can be charged only with the original charger.

The mains voltage shall conform with the specification on the charger.

Charging shall be done in accordance with section 8 of the instruction manual.

### 5 - ASPIRATING AND DISPENSING LIQUIDS

#### Attaching a pipette

**Before attaching a pipette check whether the pipette is not damaged, has no dents or sharp edges in the gripping part. Check whether the gripping part is dry.**

The pipette shall be gripped as close to the upper end as possible and carefully inserted into the pipette holder until resistance is noticed (fig. 3A). Do not apply excessive force so as not to damage thin pipettes and to avoid injury risk. A pipette that has been correctly attached and sealed in the holder should not tilt to the sides.

After attaching a pipette hold the device in such a way as to keep the pipettor in vertical position. The shape of the device casing makes it possible to put aside the pipette controller with an attached pipette. After work is finished it is recommended not to leave the device with an attached pipette for a longer period, for example overnight or over a weekend.

**Do not put aside the pipette controller if there is liquid in the pipette.**

#### Filling the pipette

Before aspirating is started, set the suction speed switch using the HIGH/LOW switch (fig. 1C):

- HIGH position - fast aspirating
- LOW position - slow aspirating

It is recommended to choose the LOW position for pipettes with a volume up to 5 ml, and the HIGH position for pipettes with a volume greater than 5 ml. Holding the pipette controller in such a way that the pipette is in vertical position immerse the pipette end in the liquid to be drawn up (fig. 3B), and press the aspiration button gently. The pipette filling speed depends on how deep the aspirating button has been pressed. The deeper the button is pressed the faster the liquid is aspirated into the pipette.

It is recommended to draw a slightly greater liquid volume than required (due to meniscus above the required volume mark), adjusting the aspiration speed, particularly in the final filling stage, so as not to overfill the pipette.

#### Setting the volume

After filling the pipette its end shall be dried with absorbent paper, that does not leave impurities, in order to remove the solution rests from the outside surface of the pipette. Then set the required liquid volume precisely. Pressing the dispense button gently (fig. 3C), dispense the excessive liquid from the pipette until the meniscus of the liquid aligns exactly with the required volume mark on the pipette.

#### Emptying the pipette

Holding the vessel in inclined position, place the pipette end in contact with the vessel wall and press the dispense button gently (fig. 3C). The dispensing intensity may be adjusted depending on how deep the dispense button has been pressed. The deeper the button is pressed the faster the outflow of liquid from the pipette.

The **Assistboy** has two dispense modes. The dispense mode is selected with the GRAV/BLOW switch (fig. 1D).

- GRAV position - the dispense is effected in gravity mode, which means that the liquid flows out of the pipette by its own weight. The outflow speed is adjusted depending on how deep the dispense button has been pressed.
- BLOW position - the dispense is effected in gravity mode like in the GRAV mode, however, when the dispense button is pressed to the end position, the pump is started and fast emptying of the pipette with a blow is effected.

#### Attention:

**During gravimetric dispensing the pipette is not completely emptied due to the characteristics of pipettes used with the pipette controller.**

### 6 - TROUBLESHOOTING

If during your work the functioning of the pipette controller is incorrect, check the cause and remove the fault.

Problem	Possible cause	Action
The pipette falls out (the holding force of the pipette is too small), or tilts to the side too much.	The pipette holder is dirty or wet (fig. 1G).	Take the pipette out of the holder. Disassemble the pipette holder and check for mechanical damage. Clean, wash and dry the pipette holder.
	The pipette holder is damaged.	If the holder is damaged, replace it with a new one.
The pump is working, but the pipette controller does not draw liquid or draws liquid very slowly.	The filter is dirty (fig. 1H).	Take the pipette out of the holder. Disassemble the pipette holder. Take out the filter - check it for impurities - if it is dirty, replace it with a new one.
	The pipette holder and/or the connector gasket is damaged (fig. 1K).	Check the pipette holder and the connector gasket for mechanical damage - if there is a damage, replace the damaged part with a new one.
Liquid leaks from the pipette (the aspiration and the dispense buttons are not pressed)	The pipette is damaged.	Check the pipette for damage - cracks, dents - if present, replace the pipette with a new one.
	The pipette is inserted incorrectly.	Check whether the pipette has been correctly inserted in the holder.
	The pipette holder, the filter, or the connector gasket is installed incorrectly.	Check whether all parts are present and correctly installed.
	The pipette holder and/or the connector gasket is damaged (fig.1G, fig. 1K).	Check the pipette holder and the connector gasket for mechanical damage - if there is damage, replace the damaged part with a new one.

If the above actions do not remedy the incorrect functioning of the pipette controller, the device shall be sent to authorised service centre.

If any manufacturing defects are discovered by the Customer, an unsatisfactory equipment claim must be compiled, certified and sent to the local distributor address. To obtain the claim form, visit section **Technical support** on our website.

## 7 - REPLACING THE FILTER

### Note:

**The work safety instructions given in section 2 shall be observed when disassembling the pipette controller.**

The filter replacement is necessary, if drawing efficiency deterioration is observed. The direct reason may be dirty filter after a long period of use.

The procedure is shown in (fig. 4).

- Remove the pipette.
- Unscrew the nosepiece (fig. 4A).
- Remove the membrane filter (fig. 4A) and the pipette holder (fig. 4B).
- Blow liquid out of the holder and leave it until it is completely dry.
- Install new membrane filter (fig. 4D) and assemble the device in reverse order.

## 8 - CHARGING THE BATTERY

**The Assistboy may be charged only with the original charger. The mains voltage must conform with the specification on the charger.**

**Using other chargers than the original one may damage the battery of the pipette controller.**

Charging the battery in the **Assistboy** pipette controller is controlled by a time circuit, which controls the entire process. When the battery has been charged, the charging circuit disconnects automatically.

### Charging:

1. Charging temperature: 10°C to 35°C.
2. Charging the battery is carried out through a charger (power supply) by direct connection to the mains, or indirectly by means of a charging stand (fig. 1M). Every **Assistboy** is supplied with a charger and a charging stand. Battery charging is indicated by successive lighting of "bars" in the display order 4-3-2-1-4-3-2-1 (fig. 5).
3. Full charging time: 11-14 hours.
4. The full battery level is indicated by displaying three "bars" which shine continuously (fig. 5.1).

The pipette controller has a NiMH type battery with a capacity of 2000 mAh.

The service life of the battery: approx. 1000 charging cycles, if used correctly. It is not possible to overcharge the battery, if all instructions of the manufacturer are followed.

### Note:

**Charging shall not be interrupted. If the charging process is interrupted, the next charging may take place after the battery has been completely discharged.**

In order to prolong the life span of a rechargeable battery the following rules should be followed:

1. Before the pipette controller is activated for the first time the battery should be charged.
2. **Assistboy** should not be left uncharged for a longer period of time.
3. If the **Assistboy** is left unused charge cycles should be performed every six months.

## 9 - MAINTENANCE

### Cleaning

The **Assistboy** does not require any maintenance. The external parts of the pipette controller may be cleaned with a swab moistened with isopropyl alcohol.

The nosepiece and the pipette holder may be autoclaved at 121°C for 20 minutes. The filters included in the set may be sterilised by autoclaving at 121°C for not more than 15 minutes.

### Ultra violet (UV) sterilization

The outer body of the pipette controller is UV resistant which was confirmed by our tests. The recommended distance from the radiation source to exposed element should be not less than 50 cm. Prolonged and very intense UV exposure can cause de-coloration of pipette controller parts, without affecting its performance.

### Storage

The **Assistboy** shall be stored in a dry place. The allowable storage temperature: -20°C to +50°C.

During breaks in the work the pipette controller shall be placed in the wall hanger included with the pipette controller.

The hanger may be fitted on a wall with a screw or stuck with the adhesive tape attached to the hanger. The surface for the hanger should be smooth, clean and degreased. After cleaning the surface, it should be wiped with isopropyl alcohol and left to dry. Before sticking the hanger remove the protective film from the tapes attached to the hanger. The hanger should be stuck by pressing it hard to the surface and by leaving it attached without any load for a period of 24 hours.

### Note:

**Do not store the pipette controller with a filled pipette.**

## 10 - COMPONENTS

The **Assistboy** set is supplied with the following components:

- Charger
- Wall hanger
- Charging stand
- Membrane filter 0.2 µm
- Instruction manual

## 11 - ORDERING INFORMATION

The **Assistboy** comes with a charger and a universal set of adapters in different versions: EU, US, UK and Australia. Choose your country used adapter and connect to the housing.

To mount the adapter, it should be inserted into the slots of the housing in the direction of the arrow (fig. 6.1), until you hear a click engagement.

To remove or change the adapter, simply press the "PUSH" button in the direction of the arrow (fig. 6.2), then holding the button down, remove the adapter in the direction of the arrow (fig. 6.3).

Catalogue no. **BS-010521**

## 12 - SPARE PARTS

Item	Part name	Cat. nr.	Number of pieces in pack
F	Nosepiece	BS-010521-S01	1
G	Pipette holder	BS-010521-S02	1
H	Membrane filter 0.2 µm	BS-010521-S03	1
	Membrane filter 0.45 µm	BS-010521-S04	1
K	Connector gasket	BS-010521-S05	1
L	Charger 9V: EU, US, UK, AU Plug	BS-010521-S06	1
M	Charging stand	BS-010521-S07	1

## 13 - PRODUCT DISPOSAL

According to Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment, pipette controller **Assistboy** is marked by the crossed-out wheeled bin. It means that the electric equipment and batteries must not be disposed of with domestic waste.

Consequently, the buyer shall follow the instructions for reuse and recycling of waste electronic and electrical equipment (WEEE).

The crossed out wheeled bin symbol is printed in the product instruction manual and on the packaging to indicate that, in accordance with the requirements of Directive 2006/66/EC of 6th September 2006 on batteries and accumulators and waste batteries and accumulators, batteries must be disposed of in accordance with national regulations.

All rights reserved. Product described in this manual is subject to availability and technical modification. Errors excepted. BIOSAN SIA reserves the right to improve, enhance or otherwise modify its products without prior notification. © 2017 BIOSAN SIA

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