

CTR-6

CO₂ tube roller



If you have any feedback on our products or services, we would like to hear from you.
Please send all feedback to:

Manufacturer:

SIA Biosan
Ratsupites 7 k-2, Riga, LV-1067, Latvia

Phone: +371 674 261 37



<https://biosan.lv>

Marketing e-mail: sales@biosan.lv

Support e-mail: support@biosan.lv

Contents

1. About this edition of the instructions	3
2. Safety precautions	4
3. General information	5
4. Getting started	6
5. Operation	10
6. Specifications	11
7. Ordering information	12
8. Care and maintenance	13
9. Storage and transportation	13
10. Warranty and registration	14
11. EU Declaration of Conformity	15

1. About this edition of the instructions

1.1. The current edition of instructions for use (IFU) applies to the following models:

Model and name	Version
CTR-6, CO ₂ tube roller	V.2AW

1.1.1. Revision number: 2.02

1.1.2. Release date: July of 2026.

1.2. Symbols used in this IFU:

 **Caution!** Please pay special attention to sections marked by this symbol.

1.3. Revision changelog

Edition	Date	Changes
2.01	05/2025	Model version change
2.02	07/2026	Logo change

2. Safety precautions

2.1. General warnings:



Caution! Make sure you have fully read and understood the present instructions before using the equipment.

2.1. Icons used on the unit and packaging:

	CE marking, manufacturer affirms conformity with European health, safety, and environmental protection standards, see 11.1
	WEEE directive marking, see 11.1
	Polarity of the power connector
	Equipment uses direct current
	Before installation in an incubator, read the user instructions!

2.2. General safety

- Save the unit from shocks or falling.
- Store and transport the unit as described in **Storage and transportation**.
- Before using any cleaning or decontamination methods except those recommended by the manufacturer, check with the manufacturer that the proposed method will not damage the equipment.
- Do not make modifications in design of the unit.
- The protection provided by the equipment may be impaired if the equipment is used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

2.3. Electrical safety

- Connect only to external power supply with voltage corresponding to that on the serial number label.
- Use only the external power supply provided with this product.
- Ensure that the external power supply is easily accessible during use.
- Disconnect the unit from the mains before moving.
- Turn off the unit by disconnecting the external power supply from the power socket.
- If liquid penetrates into the unit, disconnect it from the power socket and have it checked by a repair and maintenance technician.
- Do not operate the unit in premises where condensation can form. Operating conditions of the unit are defined in section **Specifications**.

2.4. During operation

- Hazards may occur if the equipment is used to mix flammable or explosive materials.
- Do not mix materials where the transfer of mechanical energy to glass could lead to a breakage.
- Do not use the equipment in hazardous atmospheres or with hazardous materials for which the equipment is not designed.
- Do not impede the platform motion.
- Do not operate the unit in environments with aggressive or explosive chemical mixtures. Please contact manufacturer for possible operation of the unit in specific atmospheres.
- Do not operate the unit if it is faulty or has been installed incorrectly.
- Do not use outside laboratory rooms.
- Do not place a load exceeding the maximum load value mentioned in section **Specifications**.

2.5. Biological safety

- It is the user's responsibility to carry out appropriate decontamination if hazardous material is spilt on or penetrates into the equipment.

3. General information

CO₂ Tube Roller **CTR-6** provides regulated rocking and rolling of maximum up to 6 rollers and is designed for use specifically in CO₂ incubators. CTR-6 is specifically designed for use in harsh environments such as CO₂ and humidity and provides reproducible results for cell culture growth. Possibility to remove rollers makes the unit flexible and allows for performing various procedures and techniques in various cultivation vessels. The specially designed remote controller allows for protection of electronics from CO₂ incubator environment, as well as does not interfere with the experiment.

Tube roller **CTR-6** incorporates a stepper motor with a guaranteed service life up to 10000 hours. It is possible to stack up to 3 units, saving valuable bench space. Typical applications include cells cultivation (eukaryotic, microbial) and general mixing (resuspension, viscous and liquid-solid suspensions).

4. Getting started

4.1. **Unpacking.** Remove packing materials carefully and retain them for future shipment or storage of the unit. Examine the unit carefully for any damage incurred during transit. The warranty does not cover in-transit damage. Warranty covers only the units transported in the original package.

4.2. **Complete set.** Package contents:

4.2.1. Standard set

- **CTR-6** Tube roller 1 pce.
- External power supply 1 pce.
- Remote controller 1 pce.
- Operating instructions, declaration of conformity 1 copy

4.2.2. Optional accessories, on request

- Rubber O-rings 1 set
- Stacking kit for two units 1 set
- Stacking kit for three units 1 set
- Roller 1 pce.

4.3. **Roller installation.** Install the roller (fig. 1/1) to the base mechanism of the rolling unit. Open the back fixation plate (fig. 1/2) by unscrewing both screws (fig. 1/3). Fit the roller on the mechanism shaft (fig. 1/4). The roller must be carefully aligned with other rollers horizontally. Close the fixation plate.





Figure 1. Roller installation


4.4. Stacking.

4.4.1. Stacking kit contains:

Component	Kit for 2 units	Kit for 3 units
Catalogue number	BS-010174-BK	BS-010174-CK
Connecting rods	4	8
Holder bracket	8	12
Headless hex socket screw	–	4
Hex screw	16	24
Flat washer	8	12
Lock washer	8	12
Hex key M4	1	1



Flat washer



Lock washer

4.4.2. **Acquiring a kit.** To purchase the kit, contact Biosan or your local Biosan representative, with information from the section **Ordering information**.

4.4.3. Getting started.

- Unpack the kit. Unpower and unload the units and disconnect the power cables.
- Prepare the working space. Ensure that there are no obstacles, there is enough space, the surface is clear and clean.
- Position the units so that the access to the mains plugs, the sockets and the power switches is unobstructed.

4.4.4. Stacking two units (figure 2).

- Lock the bracket parts (fig. 2/1) with hex screws and flat washers (fig. 2/2) in 4 positions on the base plate for each stacking unit.
- Insert the connecting rod (fig. 2/3) on the bottom bracket and secure it tightly with a hex screw and a lock washer. Repeat for other connecting rods.
- Place the second unit on top of the first-level rods, insert and secure with the screws and lock washers.
- Connect the power cables to the units and to grounded mains sockets.

4.4.5. Stacking three units.

- Lock the bracket parts (fig. 2/1) with hex screws and flat washers (fig. 2/2) in 4 positions on the base plate for each stacking unit.
- Insert the connecting rod (fig. 2/3) on the bottom bracket and secure it tightly with a hex screw and a lock washer. Repeat for other connecting rods.
- Screw in the headless hex socket screws into the first-level connecting rods.
- Place the second unit on top of the first unit, aligning the headless screws and brackets of the second unit.
- Screw the connecting rod onto the headless screw, securing rods and the bracket. Repeat for other connecting rods.
- Place the third unit on top of the second-level rods, insert and secure with the screws and lock washers.
- Connect the power cables to the units and to grounded mains sockets.

4.4.6. **Maintenance.** Check the screw fixation and unit stack stability once a month.




Figure 2. Stacking kit installation

4.5. Setup.



Note.

The CO₂ incubator must be not operating, switched off and in non-condensing environment during the installation of CO₂ tube roller.

- Gently pull the control cable apart at the middle connector by pressing on clips.
- Thread the control cable of the shaking unit through the opening port in the CO₂ incubator.
- Place the shaking unit inside a CO₂ incubator positioned on an even horizontal and firm surface place. Make sure that the control cable does not obstruct the movement of the platform.
- Accurately connect the control cable in the middle. Align ends by the white marking on the rim. Do not force the connection.
- Remove protective film from the display.
- Plug the external power supply into the 12 V socket at the bottom of the unit controller and position the power supply so that the plug is easily accessible.
- If necessary, install the unit controller on the outer wall or a door of the CO₂ incubator. The unit controller features magnets on the back.



Figure 3. Stacked units

5. Operation




Figure 4. Control panel of CTR-6

5.1. Connect the external power supply to the mains. Display switches on.

5.2. Place samples on the unit platform.

5.3. **Setting the parameters.** Press the **Select** key (fig. 4/1) to choose the parameter to change. Each pressing of the **Select** key consecutively activates the parameters. The active parameter is flashing. Use the ▼ and ▲ keys (fig. 4/2) to set the necessary value. Pressing the key for more than 2 s increases the speed of value change.

5.3.1. Set the required working time interval in hours and minutes, the increment is 1 minute.

5.3.2. Set the required shaking speed, the increment is 1 rpm.

5.4. Press the **Run Stop** key (fig. 4/3). The platform starts rotation, indication RUN appears on display and the timer in the lower line of the display starts counting the time interval.

5.5. After the timer reaches the set time, the platform motion will stop and the flashing indication STOP, accompanied by the repetitive sound signal, will appear in the lower line of the display. Press the **Run Stop** key to shut down the signal.

5.6. The unit can be stopped before the set time elapses, if necessary, by pressing the **Run Stop** key. Press the **Run Stop** key to repeat the operation with the same working time and speed.

5.7. If the working time is not set (or is reset) and the Time indicator on display shows OFF, pressing the **Run Stop** key will start continuous operation of the unit until the **Run Stop** key is pressed.

5.8. The platform motion can be stopped at any time by pressing the **Run Stop** key. In this case, the program realization and the platform motion will stop, and the unit will switch into the STOP mode.

5.9. **Power interruption.** In case of mains/power failure, the unit is not equipped with automatic restart function, which requires the user to reinitiate the platform motion manually.

5.10. After finishing the operation, switch the power switch on the rear side of the unit into position **O** (off) and disconnect the external power supply from the mains.

6. Specifications

6.1. Biosan is committed to a continuous programme of improvement and reserves the right to alter design and specifications of the equipment without additional notice.

6.2. Rotation specifications

Speed control range5–80 RPM



Note. Maximum speed depends on the load on the platform and the shape of the vessels

Increment 1 RPM

Digital time setting 1 min – 96 h or non-stop

Increment 1 minute

6.3. General specifications

Maximum load 3 kg

Tilt angle 4°

Dimensions, WxDxH..... 310x262x80 mm

Operating voltage and current 12 V=, 415 mA

Power consumption 5 W

External power supply input 100–240 V~, 50–60 Hz, output 12 V=

Weight, accurate within ±10% 3 kg

6.4. Workroom requirements

CTR-6 part	Shaking unit (w/o controller unit)	Controller unit
Workroom description	Specifically designed for CO ₂ incubators. Cold rooms and closed laboratory rooms, incubators	Cold rooms and closed laboratory rooms, incubators (excluding CO ₂ incubators)
Temperature range	+4 °C ... +40 °C	+4 °C ... +45 °C
Humidity requirements	Maximum of 98% RH. Non-condensing atmosphere.	Maximum of 80% RH at 31 °C, decreasing linearly to 50% RH at 40 °C. Non-condensing atmosphere.
Operating height	Maximum 2000 m ASL	
Overvoltage category	I	
Pollution degree	2	

7. Ordering information

7.1. Models and versions available:

Model	Version	Electric plug	Catalogue number
CTR-6, CO ₂ tube roller	V.2AW	230 V~, EU plug (type C)	BS-010174-A01
		100–240 V~, Multi plug	BS-010174-A02

7.2. To inquire about or order the optional accessories, contact Biosan or your local Biosan representative.

7.2.1. Optional accessories:

Description	Catalogue number
Rubber O-ring set, different diameter	BS-010174-DK
Stacking kit for two rollers	BS-010174-BK
Stacking kit for three rollers	BS-010174-CK
Roller	BS-010174-FK

8. Care and maintenance

8.1. Service.

8.1.1. If the unit is disabled (e.g., no platform movement, no reaction to key presses, etc.) or requires maintenance, disconnect the unit from the mains and contact Biosan or your local Biosan representative.

8.1.2. All maintenance and repair operations (except listed below) must be performed only by qualified and specially trained personnel.

8.1.3. Operating integrity check. If the unit follows the procedure described in section **Operation**, then no additional checks are required.

8.1.4. If the stacking kit is used, check the screw fixation monthly.

8.2. Cleaning and disinfection.

8.2.1. Use mild soap and water with a soft cloth or sponge for cleaning the exterior. Rinse remaining washing solution with distilled water. Wipe dry the excess water with clean, soft cloth or sponge.

8.2.2. To disinfect the plastic and metal parts, use 75% ethanol or DNA/RNA removing solution (e.g., Biosan PDS-250). After disinfecting, wipe the surfaces dry.

8.2.3. The unit and its accessories are not autoclavable.

8.3. **Disposal.** Disposal of the appliance requires special precautions and must be carried out at an appropriate disposal site, separate from normal household waste. To prevent pollution of the environment, all waste resulting from the disposal of the product must be collected and disposed of in the country of use, in accordance with the applicable requirements for the handling of electronic waste.

9. Storage and transportation

9.1. Store and transport the unit in a horizontal position (see package label) at ambient temperatures between -20°C and +60°C and maximum relative humidity of 80%.

9.2. Save the unit from shocks or falling.

9.3. After transportation or storage and before connecting it to the electric circuit, keep the unit under room temperature for 2-3 hrs.

9.4. For extended storage, the unit does not require special procedures.

10. Warranty and registration

- 10.1. The Manufacturer guarantees the compliance of the unit with the requirements of Specifications, provided the Customer follows the operation, storage and transportation instructions.
- 10.2. The warranted service life of the unit from the date of its delivery to the Customer is 24 months (excluding optional accessories listed in 7.2). For extended warranty, see 10.5.
- 10.3. Warranty covers only the units transported in the original package.
- 10.4. If any manufacturing defects are discovered by the Customer, an unsatisfactory equipment claim shall be compiled, certified and sent to the local distributor address. To obtain the claim form, visit section **Technical support** on our website at link below.
- 10.5. Extended warranty. For **CTR-6**, the *Basic plus* class model, extended warranty is a paid service. Please contact Biosan or your local Biosan representative.
- 10.6. Description of the classes of our products is available in the **Product class description** section on our website at the link below.

Technical support



biosan.lv/en/support

Product class description



biosan.lv/classes-en

10.7. The following information will be required in the event that warranty or post-warranty service comes necessary. Complete the table below and retain for your records.

Model	Serial number	Date of sale
CTR-6, CO ₂ tube roller		

10.8. **Production date.** Production date is placed in the serial number, on the label of the unit. Serial number consists of 14 digits styled XXXXXXYYMMZZZZ, where XXXXXX is model code, YY and MM – year and month of production, ZZZZ – unit number.

11. EU Declaration of Conformity

11.1. CO₂ tube roller **CTR-6** is in conformity with the following relevant Union legislations:

LVD 2014/35/EU	LVS EN 61010-1:2011 + A1:2019 Safety requirements for electrical equipment for measurement, control, and laboratory use. General requirements. LVS EN 61010-2-051:2015 Particular requirements for laboratory equipment for mixing and stirring.
EMC 2014/30/EU	LVS EN 61326-1:2021 Electrical equipment for measurement, control and laboratory use. EMC requirements. General requirements.
RoHS3 2015/863/EU	Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment.
WEEE 2012/19/EU	Directive on waste electrical and electronic equipment.

11.2. Declaration of Conformity is available on request.



Biosan SIA

Ratsupites 7 k-2, Riga, LV-1067, Latvia

Phone: +371 67426137

<https://biosan.lv/>

Edition 1.02 – July of 2025