Biosan SIA, Ratsupites iela 7 k-2, Riga, Latvia, LV-1067 E-mail: marketing@biosan.lv, www.biosan.lv Phone:+37167426137



Nabi, UV/VIS nano spectrophotometer

DESCRIPTION

Nabi – UV/VIS Nano Spectrophotometer can conveniently read samples in cuvette or microvolume, and also receive results from single wavelength and full spectrum Nabi is a standalone instrument which utilizes 7.0 inch touch-enabled LC. All data can be transferred to PC just with a USB.

- Sample processing in cuvette and micro volume pedestal
- Spectrometry analyzes single or full spectrum- Small and portable size
- Easy maneuver/operation with 7" LCD screen (Stand-alone)
- Easy back up with a USB Flash drive
- Parameters: DNA/RNA, Protein, Cell, End point, Kinetic, Spectrum

More info



Nabi, UV/VIS nano spectrophotometer Page 1 of 1