




Use only tubes, plates suitable for centrifugation, adhere to the vessel manufacturer's specifications for optimal centrifugation parameters, including load and speed.




## STANDARD ROUND BOTTOM CENTRIFUGE TUBE



**10-15 ml**  
Ø: 16 mm  
L: 105 mm


- Ri-12/10**  
4,200 rpm  
BS-010221-MK
- Ri-24/10**  
4,000 rpm  
BS-010221-LK
- BR-4U + BA-4/RT**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-3K <sup>+</sup>

## VACUTAINERS




**2-5 ml**  
Ø: 13 mm  
L: 82 mm  
[13 × 75]

- Ri-12/10**  
4,200 rpm  
BS-010221-MK
- Ri-24/10**  
4,000 rpm  
BS-010221-LK
- BR-4U + BA-4/VT**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-4K <sup>+</sup>



**4-8 ml**  
Ø: 13 mm  
L: 107 mm  
[13 × 100]


- Ri-12/10**  
4,200 rpm  
BS-010221-MK
- Ri-24/10**  
4,000 rpm  
BS-010221-LK
- BR-4U + BA-4/VT**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-4K <sup>+</sup>



**8-9 ml**  
Ø: 16 mm  
L: 105 mm  
[16 × 100]


- Ri-12/10**  
4,200 rpm  
BS-010221-MK
- Ri-24/10**  
4,000 rpm  
BS-010221-LK
- BR-4U + BA-4/VT**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-4K <sup>+</sup>

## STANDARD CONICAL BOTTOM CENTRIFUGE TUBE




**15 ml**  
Ø: 17 mm  
L: 120 mm

- Ri-12/15**  
4,200 rpm  
BS-010221-KK
- BR-4U + BA-4/FT**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-YK <sup>+</sup>



**25 ml**  
Ø: 29 mm  
L: 77 mm

- Ri-6**  
4,200 rpm  
BS-010221-HK
- Ri-6P**  
4,200 rpm  
BS-010221-IK
- BR-4U + BA-2/50**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-CK <sup>+</sup>



**50 ml**  
Ø: 29 mm  
L: 115 mm

- Ri-6**  
4,200 rpm  
BS-010221-HK
- Ri-6P**  
4,200 rpm  
BS-010221-IK
- BR-4U + BA-2/50**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-CK <sup>+</sup>

## MICROTUBES



**1.5-2 ml**  
Ø: 11 mm  
L: 40 mm

- RMT-24**  
6,000 rpm  
BS-010221-BK



**5 ml**  
Ø: 16 mm  
L: 60 mm


- BR-4U + BA-14/2U**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-JK <sup>+</sup>
- BR-4U + BA-8/5**  
4,200 rpm  
BS-010221-AK <sup>+</sup>  
BS-010221-FK <sup>+</sup>



Use only tubes, plates suitable for centrifugation, adhere to the vessel manufacturer's specifications for optimal centrifugation parameters, including load and speed.



## 96, 384-WELL, DEEP-WELL PLATES, MICROPLATES



→ (x2) →




**Ri-2**  
2,000 rpm  
BS-010221-GK


Plate dimension

L: 128 mm W: 85.6 mm H<sub>max</sub>: 45mm

## SEMI- AND UNSKIRTED 96-WELL MICROPLATES



→ (x2) →




**Ri-2**  
2,000 rpm  
BS-010221-GK

Plate dimension

L: 128 mm W: 85.6 mm H<sub>max</sub>: 45mm

## UNSKIRTED 384-WELL MICROPLATES




**Ri-2**  
2,000 rpm  
BS-010221-GK

Plate dimension

L: 128 mm W: 85.6 mm H<sub>max</sub>: 45mm

## GEL CARDS



→ (x24) →



**Ri-24GC**  
1,500 rpm  
BS-010221-NK

Dimension

L: 74 mm H: 53 mm

## RACK FOR ROTORS



→ (x1) →



**RR-U**  
Rack for rotors  
BS-010208-UK