



50ml - PSU-20i - 42 pcs




UP-168

50ml - ES-20/60, ES-20/80, ES-20/80C - 36 pcs




UP-168

100ml - PSU-20i - 20 pcs




UP-168

100ml - ES-20/60, ES-20/80, ES-20/80C - 16 pcs




UP-168

250ml - PSU-20i - 14 pcs




UP-168

250ml - ES-20/60, ES-20/80, ES-20/80C - 11pcs




UP-168

500ml - PSU-20i - 12 pcs




UP-168

500ml - ES-20/60, ES-20/80, ES-20/80C - 9pcs




UP-168

1000ml - PSU-20i - 8 pcs




UP-168

1000ml - ES-20/60, ES-20/80, ES-20/80C - 5pcs




UP-168

2000ml - PSU-20i - 4 pcs




UP-168

2000ml - ES-20/60, ES-20/80, ES-20/80C - 4 pcs
Using P-EX, platform exchange set.




UP-168

TR-44/15, TR-21/50 - PSU-20i, ES-20/60, ES-20/80, ES-20/80C - 2 pcs



UP-168

SPML - PSU-20i, ES-20/60, ES-20/80, ES-20/80C - 4 pcs



UP-168