

# UP-12, platform

## DESCRIPTION

Universal platform with adjustable bars for different types of flasks, bottles, and beakers with silicone mat.

## SPECIFICATIONS

Working area	270 × 185 × 40 mm
Flask volume	universal



## CAT. NUMBER

BS-010108-AK

Platform UP-12

## ACCESSORIES



**HB-200**  
 BS-010108-FK

Additional holding bar for UP-12



**ES-20**  
 Without platform  
 Orbital Shaker-Incubator

The **ES-20** is a compact benchtop Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



**PSU-10i**  
 Without platform  
 Orbital Shaker

The **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



**CPS-20**  
 Without platform  
 CO<sub>2</sub> Shaker

The **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in CO<sub>2</sub> incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)