

P-12/100, platform

DESCRIPTION

Platform with clamps for flasks.

SPECIFICATIONS

Working area	250 × 190 mm
Flask volume	100-150 ml
Number of clamps	12



CAT. NUMBER

BS-010108-EK

Platform P-12/100

ACCESSORIES



ES-20
 Without platform
 Orbital Shaker-Incubator

The **ES-20** is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



PSU-10i
 Without platform
 Orbital Shaker

The **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



CPS-20
 Without platform
 CO₂ Shaker

The **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in CO₂ incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)