


## HSC-50, clamp

### DESCRIPTION

Tight fit clamp 50 ml - Ø50 mm for UP-168



### CAT. NUMBER

BS-010167-DK

HSC-50, clamp 50ml

### ACCESSORIES



**UP-168**  
BS-010135-JK  
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



**ES-20/80**  
Software included, without  
platform  
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



**ES-20/80C**  
Software included, without  
platform  
Orbital Shaker-Incubator with  
cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)