

P-05/32, platform

DESCRIPTION

For 32 tubes 0.5 ml.

SPECIFICATIONS

Capacity	32
Microtest tubes volume	0.5 ml



CAT. NUMBER

BS-010216-BK	Platform P-05/32
--------------	------------------



MPS-1

With built-in universal platform
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)