

# PRS-1DP, platform

## DESCRIPTION

Platform for microplates and racks for tall tubes 0.5 and 1 ml (e.g. Thermo 3741MTX, 3742MTX, 3744MTX).

## SPECIFICATIONS

Tube volume

0.5 ml and 1 ml



## CAT. NUMBER

BS-010149-DK

Platform PRS-1DP



**Multi Bio RS-24**  
Including PRS-26 platform  
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)