

# SC-18, block

## DESCRIPTION

20 × 0.5 ml + 12 × 1.5 ml microtubes

## SPECIFICATIONS

Amount of 0.5 ml tubes	20
Amount of 1.5 ml tubes	12



## CAT. NUMBER

BS-010120-AK	SC-18
--------------	-------



#### TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)