

## SC-24C, block

### DESCRIPTION

24 × 2 ml microtubes

### SPECIFICATIONS

Amount of 2 ml tubes

24



### CAT. NUMBER

BS-010143-EK

SC-24C

### ACCESSORIES



#### TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



#### TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)