

MSL-SC, protection lid

DESCRIPTION

Higher rotor lid to accommodate spin columns



CAT. NUMBER

BS-010213-EK

MSL-SC, higher rotor lid to accommodate spin columns



Microspin 12

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical analysis of ...

[read more](#)