


MA-8, 8-channel adapter manifold

DESCRIPTION

8-channel aspiration tip with holder, autoclavable



CAT. NUMBER

BS-040108-BK

MA-8, Optional 8-channel adapter manifold, MA-8, 8-channel adapter manifold

ACCESSORIES



FTA-1
BS-040108-AAG
Aspirator with trap flask

Aspirator with trap flask **FTA-1** is designed for aspiration/removal of alcohol/buffer remaining quantities from microtest tube walls during DNA/RNA purification and other macromolecule reprecipitation techniques.

The device can be used ...

[read more](#)



FTA-2i
BS-040120-A02
Aspirator with Trap Flask

Aspirator with trap flask **FTA-2i** is designed for aspiration or removal of alcohol, buffer and liquid from reaction vessels (e.g. during DNA/RNA purification or other macromolecule reprecipitation techniques).

The device ...

[read more](#)