

1-position stand

DESCRIPTION

SPECIFICATIONS

Capacity

1



CAT. NUMBER

BS-010524

1-position stand

ACCESSORIES



Assist

Adjustable volume pipettes - Single channel
 pipette series

The **Assist series** pipettes are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on the model.

Multichannel ...

[read more](#)



Assist

Adjustable volume pipettes - Single channel
 pipette series with volume lock

The **Assist series** pipettes with volume lock are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on

...

[read more](#)