

IPP-4, platform

DESCRIPTION

Platform IPP-4


SPECIFICATIONS

Working area

266 x 170 mm

Flask volume

4 microtest plates



CAT. NUMBER

BS-010102-AK

Platform IPP-4



PSU-2T
Including IPP-2 platform
Mini-shaker for immunology

Mini-Shaker **PSU-2T** is designed for immunoassays and provides adjustable mixing of reagents in microplates. The device ensures smooth movement of the platform even at low speeds.

Shaker is a compact ...

[read more](#)