


Extended tubing

DESCRIPTION

Extended tubing 2 m long, with fittings and MA-U adapter



CAT. NUMBER

BS-040120-DK

Extended tubing 2 m long, with fittings and MA-U adapter



FTA-2i
BS-040120-A02
Aspirator with Trap Flask

Aspirator with trap flask **FTA-2i** is designed for aspiration or removal of alcohol, buffer and liquid from reaction vessels (e.g. during DNA/RNA purification or other macromolecule reprecipitation techniques).

The device ...

[read more](#)