

1000 µl filter tips

DESCRIPTION

Pipette tips features:

- The tip is made of high-quality medical grade polypropylene, and the filter element is made of ultra-high molecular weight polyethylene, which has good hydrophobicity
- The filter element effectively protects the pipette from the sample to ensure the safety of sample aspiration
- No inner surface coating, will not pollute the sample
- Have precise tick marks
- DNase, RNase and Pyrogen - free
- Suitable for conventional pipettes such as Biosan (Assist series) Eppendorf and Gilson
- Provided in sterile racks, 96 pcs per rack



SPECIFICATIONS

Dimensions (H×Ø)

86 × 0.8 mm

CAT. NUMBER

BHZ05R1W-FS

1000 µl - 96 pcs in stand

ACCESSORIES



Assist

Adjustable volume pipettes - Single channel
pipette series

The **Assist series** pipettes are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on the model.

Multichannel ...

[read more](#)



Assist

Adjustable volume pipettes - Single channel
pipette series with volume lock

The **Assist series** pipettes with volume lock are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on

...

[read more](#)