

Gentier 96E, Real-time PCR system

DESCRIPTION

The Tianlong Gentier 96E real-time PCR system is a PCR test machine designed to meet the specific needs of high-end laboratories. Featuring six fluorescence channels, the Gentier 96E is capable of processing 96 samples in one test run. Indulge in the powerful and efficient temperature control system, easy-to-use software, and user-friendly operational designs. Tianlong Gentier 96E real-time PCR system can provide maximal reliability and efficiency for all your real-time PCR needs.

[Read more](#)



CAT. NUMBER

Gentier 96E

230VAC 50/60Hz Euro plug