

## TSP-24/1.5, block

### BESCHREIBUNG

24 x 1,5ml-Mikroröhrchen

### KAT.-NUMMER

BS-010178-KK

TSP-24/1.5

### TECHNISCHE DATEN

|                           |       |
|---------------------------|-------|
| Menge von 1,5 ml Röhrchen | 24    |
| Maximale Drehzahl, U/min  | 2500  |
| Maximaltemperatur         | 100°C |

### ZUBEHÖR



#### TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)