

TSP-24/2.0, Block

BESCHREIBUNG

24 × 2 ml-Mikroröhrchen

KAT.-NUMMER

BS-010178-JK

TSP-24/2.0

TECHNISCHE DATEN

Menge von 2 ml Röhrchen	24
Maximale Drehzahl, U/min	2500
Maximaltemperatur	100°C

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)