

TSP-8/5, Block

BESCHREIBUNG

8 × 5 ml Röhrchen

KAT.-NUMMER

BS-010178-IK

TSP-8/5

TECHNISCHE DATEN

Maximale Drehzahl, U/min

1200

Maximaltemperatur

100°C

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)