

TSP-4/50, block

BESCHREIBUNG

4 × 50 ml Röhrchen

KAT.-NUMMER

BS-010178-GK

TSP-4/50

TECHNISCHE DATEN

Maximale Drehzahl, U/min	1000
Maximaltemperatur	80°C
Ø	29.5 mm
Tiefe	76.5 mm

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)