

TSP-20/CV10, Block

BESCHREIBUNG

20 × 10mm Küvetten

KAT.-NUMMER

BS-010178-DK

TSP-20/CV10

TECHNISCHE DATEN

Maximale Drehzahl, U/min	1100
--------------------------	------

Maximaltemperatur	80°C
-------------------	------

∅	□ 12.5 x 12.5 mm
---	------------------

Tiefe	31 mm
-------	-------

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)