

TSP-24/CRYO4, Block

BESCHREIBUNG

24 × 3,6–4,5 ml Kryoröhrchen

KAT.-NUMMER

BS-010178-CK

TSP-24/CRYO4

TECHNISCHE DATEN

Maximale Drehzahl, U/min	1300
Maximaltemperatur	100°C
Ø	12.5 mm
Tiefe	30 mm

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The TS-100C Pro is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)