

TSP-24/CRY01.8, Block

BESCHREIBUNG

24 × 1–1,8 ml Kryoröhrchen

KAT.-NUMMER

BS-010178-BK

TSP-24/CRY01.8

TECHNISCHE DATEN

Maximale Drehzahl, U/min	1400
Maximaltemperatur	100°C
Ø	12.5 mm
Tiefe	20 mm

ZUBEHÖR



TS-100C Pro

Programmierbarer Thermoschüttler mit Heiz- und Kühlfunktion

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Mehr](#)