

Sistema para preparación de medios SMP-160

DESCRIPCIÓN

System of medium Preparation **SMP-160** designed for Levinsteina-Jensena medium preparation for the cultivation of mycobacterium. Instrument may be used for the heating of various reagents too.

Instrument has digital setting and indication of the temperature and time interval of incubation in stand.

Available two versions with tube holder or with shelves.

ESPECIFICACIONES

Intervalo de regulación de temperatura	+30..+95°C
Estabilidad de temperatura	±0.1 °C
Ajuste de temperatura	digital
Ajuste de tiempo digital	1 min - 96 hrs
Volumen de cámara	45 l
Número de sitios en el soporte	186
Dimensión del tubo (longitud x diámetro)	110 mm x 16 mm
Desviación de la temperatura de la cámara con respecto a la temperatura fijada	no more than 2%
Dimensiones generales (An. x Prof. x Al.)	400x410x580 mm
Peso	51 kg
Alimentación consumida, no superior a	600 W
Tensión nominal	230V; 50/60Hz, 115V; 50/60Hz



NÚMERO DE CATEGORÍA

BS-010136-A01

Enchufe europeo de 230 V CA 50/60 Hz | con estantes

BS-010136-BAA

230VAC 50/60Hz Euro plug | with tube rack