

TSP-24/2.0, bloc

DESCRIPTION

24 microtubes de 2 ml

CARACTÉRISTIQUES

Quantité de tubes de 2 ml	24
RPM maximum, tr/min	2500
Température maximale	100°C

NUMÉRO DE CAT.

BS-010178-JK

TSP-24/2.0

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)