

TSP-8/5, bloc

DESCRIPTION

8 tubes coniques de 5 ml

CARACTÉRISTIQUES

RPM maximum, tr/min

1200

Température maximale

100°C

NUMÉRO DE CAT.

BS-010178-IK

TSP-8/5

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)