

TSP-8/15, bloc

DESCRIPTION

8 tubes coniques de 15 ml

CARACTÉRISTIQUES

RPM maximum, tr/min	1100
Température maximale	80°C
Ø	17 mm
Profondeur	80 mm

NUMÉRO DE CAT.

BS-010178-FK

TSP-8/15

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The TS-100C Pro is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)