

TSP-4/50, bloc

DESCRIPTION

4 tubes coniques de 50 ml

CARACTÉRISTIQUES

RPM maximum, tr/min	1000
Température maximale	80°C
Ø	29.5 mm
Profondeur	76.5 mm

NUMÉRO DE CAT.

BS-010178-GK

TSP-4/50

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The TS-100C Pro is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)