

TSP-20/CV10, bloc

DESCRIPTION

20 cuvettes de 10 mm

CARACTÉRISTIQUES

RPM maximum, tr/min	1100
Température maximale	80°C
Ø	□ 12.5 x 12.5 mm
Profondeur	31 mm

NUMÉRO DE CAT.

BS-010178-DK

TSP-20/CV10

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)