

TSP-24/CRY01.8, bloc

DESCRIPTION

24 cryotubes de 1 à 1,8 ml

CARACTÉRISTIQUES

RPM maximum, tr/min	1400
Température maximale	100°C
Ø	12.5 mm
Profondeur	20 mm

NUMÉRO DE CAT.

BS-010178-BK

TSP-24/CRY01.8

ACCESSOIRES



TS-100C Pro

Thermo-agitateur programmable avec chauffage et refroidissement

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lire la suite](#)