

TSP-24/CRY01.8, blocco

DESCRIZIONE

24 criotubi da 1–1,8 ml

SPECIFICHE

RPM massimo	1400
Temperatura massima	100°C
Ø	12.5 mm
Profondità	20 mm

NUMERO DI CATEGORIA

BS-010178-BK

TSP-24/CRY01.8

ACCESSORI



TS-100C Pro

Termoagitatore programmabile con riscaldamento e raffreddamento

The TS-100C Pro is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Leggi di più](#)