

# qPCR (PKR) reaģenti

## APRAKSTS

Produkti:

- [BlooDirect TaqProbe qPCR Master Mix, Low ROX, 4x](#)
- [CitriNStart Taq PCR Master Mix, 2x](#)
- [EmeraldNStart HiFi Marathon PCR Master Mix](#)
- [First-Strand cDNA Synthesis Master Mix, 4x](#)
- [GoldNStart TaqGreen qPCR Master Mix](#)
- [Lyophilized GoldNStart TaqGreen qPCR Master Mix](#)
- [Nuclease-free Water, 20 mL](#)
- [One-Step GoldNStart TaqGreen RT-qPCR Master Mix, 4x](#)
- [One-Step TaqProbe RT-qPCR Master Mix, 4x](#)
- [qPCR Loading Indicator Kit, 1,500 reactions](#)
- [SapphiNStart TaqProbe qPCR Master Mix](#)
- [Single Cell Whole Genome Amplification Kit](#)

