

TSP-24/CRY04, bloks

APRAKSTS

24 × 3.6–4.5 ml kriomēģenēm

KAT. NUMURS

BS-010178-CK

TSP-24/CRY04

SPECIFIKĀCIJA

Maksimālie apgr./min	1300
Maksimālā temperatūra	100°C
Bloka ligzdas $\varnothing$	12.5 mm
Bloka ligzdas dziļums	30 mm

AKSESUĀRI



TS-100C Pro

Programmējams termokratītājs ar sildīšanu un dzesēšanu

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lasīt vairāk](#)