

TSP-24/CRY01.8, bloks

APRAKSTS

24 × 1–1.8 ml kriomēģenēm

SPECIFIKĀCIJA

Maksimālie apgr./min	1400
Maksimālā temperatūra	100°C
Bloka ligzdas $\varnothing$	12.5 mm
Bloka ligzdas dziļums	20 mm

AKSESUĀRI

KAT. NUMURS

BS-010178-BK

TSP-24/CRY01.8



TS-100C Pro

Programmējams termokratītājs ar sildīšanu un dzesēšanu

The TS-100C Pro is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[Lasīt vairāk](#)