

T-4, Table

DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

SPECIFICATIONS

Maximum load	50 kg
Overall dimensions (W×D×H)	800 × 600 × 770 mm
Weight	23 kg



CAT. NUMBER

BS-040101-BK

T-4, Table

ACCESSORIES



UVC/T-AR
BS-040102-A13
PCR cabinet

The Biosan **UVC/T-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVC/T-M-AR
With built in socket and inlet for power cords
PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVT-B-AR
With built in socket and inlet for power cords
PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



LF-1
BS-050101-BK
Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.