

# P-12/100, platform

## DESCRIPTION

Platform with clamps for flasks.

## SPECIFICATIONS

Working area	250 × 190 mm
Flask volume	100-150 ml
Number of clamps	12



## CAT. NUMBER

BS-010108-EK

Platform P-12/100

## ACCESSORIES



ES-20

Without platform  
Orbital Shaker-Incubator

The ES-20 is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



PSU-10i

Without platform  
Orbital Shaker

Shaker PSU-10i provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of six (6)

...

[read more](#)



CPS-20

Without platform  
CO<sub>2</sub> Shaker

CO<sub>2</sub> Shaker CPS-20 provides regulated orbital motion of the platform and is designed for use specifically in CO<sub>2</sub> incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)