

# P-16/88, platform

## DESCRIPTION

Platform with spring holders for up to 88 tubes up to 30 mm diameter (e. g. 10 ml, 15 ml, 50 ml tubes).

## SPECIFICATIONS

Working area	275 × 205 × 75 mm
--------------	-------------------



## CAT. NUMBER

BS-010116-BK	Platform P-16/88
--------------	------------------

## ACCESSORIES



**ES-20**  
 Without platform  
 Orbital Shaker-Incubator

The **ES-20** is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



**PSU-10i**  
 Without platform  
 Orbital Shaker

The **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



**CPS-20**  
 Without platform  
 CO<sub>2</sub> Shaker

The **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in CO<sub>2</sub> incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)