

P-16/250, platform

DESCRIPTION

Platform with clamps for flasks.

SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	250-300 ml
Number of clamps	16



CAT. NUMBER

BS-010135-CK

Platform P-16/250

ACCESSORIES



ES-20/60
Without platform
Orbital Shaker-Incubator

Orbital Shaker-Incubator ES-20/60 for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



PSU-20i
Without platform
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)