

HSP-6/1000, platform

DESCRIPTION

Platform with tight fit clamps for flasks.

SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	1000 ml
Number of clamps	6




CAT. NUMBER

BS-010167-LK

Platform HSP-6/1000

ACCESSORIES



P-EX
BS-010173-CK
platform exchange set

Platform Exchange Set allows for faster replacement of flat platforms. Thanks to the platform sliding system, it is now possible to install 4 × 2l flasks in ES-20/60, ES-20/80, ES-20/80C ...

[read more](#)



ES-20/80
Software included, without platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C
Software included, without platform
Orbital Shaker-Incubator with cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)