

FC-1000, clamp

DESCRIPTION

Clamp 1000 ml - Ø130 mm for UP-168



CAT. NUMBER

BS-010126-IK

FC-1000, clamp 1000ml

ACCESSORIES



UP-168
BS-010135-IK
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



ES-20/60
Without platform
Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



PSU-20i
Without platform
Multi-functional Orbital Shaker

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)