

# FC-1000, clamp

## DESCRIPTION

Clamp 1000 ml - Ø130 mm for UP-168



## CAT. NUMBER

BS-010126-IK

FC-1000, clamp 1000ml

## ACCESSORIES



**UP-168**  
BS-010135-JK  
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



**ES-20/60**  
Without platform  
Orbital Shaker-Incubator

Orbital Shaker-Incubator **ES-20/60** for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



**PSU-20i**  
Without platform  
Multifunctional Orbital Shaker

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)