

HSC-500, clamp

DESCRIPTION

Tight fit clamp 500 ml - Ø105 mm for UP-168



CAT. NUMBER

BS-010167-JK

HSC-500, clamp 500ml

ACCESSORIES



UP-168
BS-010135-JK
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



ES-20/80
Software included, without
platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C
Software included, without
platform
Orbital Shaker-Incubator with
cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)