

P-2/24, platform

DESCRIPTION

For 24 tubes 1.5-2 ml.

SPECIFICATIONS

Capacity	24
Microtubes volume	1.5-2 ml



CAT. NUMBER

BS-010216-AK Platform P-2/24



MPS-1
With built-in universal platform
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)