

PRS-26, platform

DESCRIPTION

Can accommodate 26 tubes 10-16 mm diameter (1.5-15 ml).

SPECIFICATIONS

Platform type	standard
Tube diameter	10-16 mm
Tube volume	1.5-15 ml
Capacity	26



CAT. NUMBER

BS-010117-GK

Platform PRS-26



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)