

PRS-5/12, platform

DESCRIPTION

5 tubes 20-30 mm diameter (50 ml tubes) and 12 tubes 10-16 mm diameter (1.5-15 ml tubes).

SPECIFICATIONS

Platform type	optional
Tube diameter	20-30 and 10-16 mm
Tube volume	up to 50 and 1.5-15 ml
Capacity	5 and 12



CAT. NUMBER

BS-010117-HK

Platform PRS-5/12



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)