

PRS-10, platform

DESCRIPTION

10 tubes 20-30 mm diameter (up to 50 ml tubes).

SPECIFICATIONS

Platform type	standard
Tube diameter	20-30 mm
Tube volume	up to 50 ml
Capacity	10



CAT. NUMBER

BS-010117-IK

Platform PRS-10



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)