

PRSC-22, platform

DESCRIPTION

Can accommodate 22 tubes 16 mm diameter (15 ml).

SPECIFICATIONS

Capacity	22
Platform type	optional
Tube diameter	25-30 mm
Tube volume	15 ml



CAT. NUMBER

BS-010117-LK

Platform PRSC-22



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)